

UNIVERSAL BOND II

The bedrock of dental restorations – standardised application on any surface.



No light curing


UNIVERSAL BOND II reacts by contact curing (BoSE technology). This eliminates the usual light curing and saves time during application.

Standardised work

No matter which surface – just 3 steps apply: Mix – Apply – Air blow – DONE! This minimises errors.

Visual control

The cleverly coloured liquids show the user the individual steps via colour change – full control.



Full visual control of the individual work steps thanks to the unique color change



The new UNIVERSAL BOND II from Tokuyama Dental can rightly be called universal bonding. It is easy to handle and achieves comparable results to previous bonding systems with fewer working steps. In combination with ESTEC EM II PLUS, it was possible to realise an aesthetically good restoration for our patient.”

Gordon Kautzsch, Leipzig