




Our RA660 range offers up-to-date radiation protection clothing which accompanies your active working day and provides outstanding service. The extraordinarily good wearing comfort is supported by the wide stretch insert on the back panels. Our customers report that this model is their favourite front apron. This is because it can be taken on and off easily, fits the body well, stays in place well and provides shoulder relief.

**Size choice:** The size depends on the girth of the chest, waist, or hip, depending on which of the three measurements is the largest.

**Note:** *Shoulder pads are also included with this RA660 model.*

**Order Details / Information**

Example for the formation of the article code:



The listed article codes are standard configurations. If you choose different options, simply change the code at the respective position.

COMFORT  
RA660

Lead Equivalent – Front	► Code	Lead Equivalent – Back	► Code
Total value or for the Levio Edition base layer / total value:		Back side:	
<b>0.25 mm Pb</b>	► 2	<b>0.25 mm Pb</b> (standard)	► 2
<b>0.35 mm Pb</b>	► 3	<b>0.00 mm Pb</b> (optional)	► 0
<b>0.50 mm Pb</b>	► 5		
<b>0.25 mm Pb / 0.42 mm Pb (Levio)</b>	► 6	Discount for no protection / 0,00 mm Pb	- 10.00 €
<b>0.25 mm Pb / 0.50 mm Pb (Levio)</b>	► 7		

Sizes, recommended for following girths:

Size	Breast	Waist	Hip	► Code
<b>XS</b>	Breast ≤ 84 cm	Waist ≤ 91 cm	Hip ≤ 91 cm	► Z
<b>S</b>	Breast ca. 85 - 95 cm	Waist ca. 92 - 102 cm	Hip ca. 92 - 102 cm	► S
<b>M</b>	Breast ca. 96 - 106 cm	Waist ca. 103 - 112 cm	Hip ca. 103 - 112 cm	► M
<b>L</b>	Breast ca. 107 - 117 cm	Waist ca. 113 - 123 cm	Hip ca. 113 - 123 cm	► L
<b>XL</b>	Breast ca. 118 - 128 cm	Waist ca. 124 - 134 cm	Hip ca. 124 - 134 cm	► X
<b>XXL</b>	Breast ca. 129 - 139 cm	Waist ca. 135 - 145 cm	Hip ca. 135 - 145 cm	► D
<b>3XL</b>	Breast ca. 140 - 150 cm	Waist ca. 146 - 156 cm	Hip ca. 146 - 156 cm	► T

Lengths, recommended height:	► Code
<b>100 cm</b> for ca. 160 - 170 cm	► 10
<b>110 cm</b> for ca. 170 - 180 cm	► 11
<b>120 cm</b> for ca. 175 - 190 cm	► 12
<b>130 cm</b> for ca. 185 - 195 cm	► 13
<b>90 cm*</b> for ca. 150 - 160 cm (Price 90 cm = Price 100 cm - 20.00 €)	► 90
<b>140 cm*</b> for ca. 195 - 205 cm (Price 140 cm = Price 130 cm + 20.00 €)	► 14

\* Special lengths, on request

Colours	► Code
<b>Curacao</b>	► 22
<b>Lagoon</b>	► 25
<b>Orchidee</b>	► 27
<b>Regatta</b>	► 28
<b>Titan</b>	► 29
<b>Mandarin</b>	► 30
<b>Blue Admiral</b>	► 31
<b>Baccara</b>	► 32

Outer Material – Cover	► Code
<b>ComforTex HPMF</b>	► 1
<b>ComforTex HPMF Hybrid</b>	► 2
Discount for ComforTex HPMF (apron)	- 20.00 €

X-Ray Protective Material – Inlay, for lead equivalents of:	► Code
<b>Standard Lead</b> 0.25 mm Pb, 0.35 mm Pb, 0.42 mm Pb and 0.50 mm Pb	► 2
<b>NovaLite</b> 0.25 mm Pb, 0.42 mm Pb and 0.50 mm Pb	► 1
<b>NovaLite</b> 0.35 mm Pb	► 8
<b>Leadfree</b> 0.25 mm Pb, 0.35 mm Pb, 0.42 mm Pb and 0.50 mm Pb	► 3

**Note: 0.42 mm Pb only for Levio Edition!**

Optional	Price/€	Additions/Replacement	Price/€
<b>Outlast – Cooling material, optional</b>	30.00	<b>Bolero</b> (see page 36)	from 185.00
<b>Sleeves matching apron</b> (see page 37)	from 85.00	<b>Thyroid and Sternum Shield</b> (see page 38)	from 94.00
<b>Breast pocket without embroidery</b>	30.00		
<b>Breast pocket with embroidery</b>	50.00		
<b>Permanent embroidery</b> (text max. three lines)	20.00		
<b>Removable embroidery</b> (text max. three lines)	20.00		
<b>Larger Sizes</b>			
• Size XXL (apron)	Price XL + 10.00		
• Size 3XL (apron)	Price XL + 20.00		
<b>Customization</b>			
• Individual (personal measurements)	On request		

**COMFORT RA660** Standard Collection with different **X-Ray Protective Materials**, a lead equivalent in the **Back of 0.25 mm Pb** (Apron) and **ComforTex HPMF Hybrid** outer material

**RA**

	Front (Pb)	Length	Article Code (Start)	Colour Code (CC)	Code Standard Lead (50 - 150 kV)	Price/€	Code NovaLite (50 - 110 kV)	Price/€	Code Leadfree (50 - 110 kV)	Price/€
Size XS	0.25 mm	100 cm	RA66022Z10	+CC	+22	310.00	+21	400.00	+23	370.00
	0.25 mm	110 cm	RA66022Z11	+CC	+22	330.00	+21	420.00	+23	390.00
	0.25 mm	120 cm	RA66022Z12	+CC	+22	350.00	+21	440.00	+23	410.00
	0.25 mm	130 cm	RA66022Z13	+CC	+22	370.00	+21	460.00	+23	430.00
	0.35 mm	100 cm	RA66032Z10	+CC	+22	370.00	+28	440.00	+23	420.00
	0.35 mm	110 cm	RA66032Z11	+CC	+22	390.00	+28	460.00	+23	440.00
	0.35 mm	120 cm	RA66032Z12	+CC	+22	410.00	+28	480.00	+23	460.00
	0.35 mm	130 cm	RA66032Z13	+CC	+22	430.00	+28	500.00	+23	480.00
Size S	0.50 mm	100 cm	RA66052Z10	+CC	+22	420.00	+21	580.00	+23	520.00
	0.50 mm	110 cm	RA66052Z11	+CC	+22	440.00	+21	600.00	+23	540.00
	0.50 mm	120 cm	RA66052Z12	+CC	+22	460.00	+21	620.00	+23	560.00
	0.50 mm	130 cm	RA66052Z13	+CC	+22	480.00	+21	640.00	+23	580.00
Size M	0.25 mm	100 cm	RA66022S10	+CC	+22	320.00	+21	410.00	+23	380.00
	0.25 mm	110 cm	RA66022S11	+CC	+22	340.00	+21	430.00	+23	400.00
	0.25 mm	120 cm	RA66022S12	+CC	+22	360.00	+21	450.00	+23	420.00
	0.25 mm	130 cm	RA66022S13	+CC	+22	380.00	+21	470.00	+23	440.00
	0.35 mm	100 cm	RA66032S10	+CC	+22	380.00	+28	450.00	+23	430.00
	0.35 mm	110 cm	RA66032S11	+CC	+22	400.00	+28	470.00	+23	450.00
	0.35 mm	120 cm	RA66032S12	+CC	+22	420.00	+28	490.00	+23	470.00
	0.35 mm	130 cm	RA66032S13	+CC	+22	440.00	+28	510.00	+23	490.00
Size L	0.50 mm	100 cm	RA66052S10	+CC	+22	430.00	+21	590.00	+23	530.00
	0.50 mm	110 cm	RA66052S11	+CC	+22	450.00	+21	610.00	+23	550.00
	0.50 mm	120 cm	RA66052S12	+CC	+22	470.00	+21	630.00	+23	570.00
	0.50 mm	130 cm	RA66052S13	+CC	+22	490.00	+21	650.00	+23	590.00
	0.25 mm	100 cm	RA66022M10	+CC	+22	330.00	+21	420.00	+23	390.00
	0.25 mm	110 cm	RA66022M11	+CC	+22	350.00	+21	440.00	+23	410.00
	0.25 mm	120 cm	RA66022M12	+CC	+22	370.00	+21	460.00	+23	430.00
	0.25 mm	130 cm	RA66022M13	+CC	+22	390.00	+21	480.00	+23	450.00
Size XL	0.35 mm	100 cm	RA66032M10	+CC	+22	390.00	+28	460.00	+23	440.00
	0.35 mm	110 cm	RA66032M11	+CC	+22	410.00	+28	480.00	+23	460.00
	0.35 mm	120 cm	RA66032M12	+CC	+22	430.00	+28	500.00	+23	480.00
	0.35 mm	130 cm	RA66032M13	+CC	+22	450.00	+28	520.00	+23	500.00
	0.50 mm	100 cm	RA66052M10	+CC	+22	440.00	+21	600.00	+23	540.00
	0.50 mm	110 cm	RA66052M11	+CC	+22	460.00	+21	620.00	+23	560.00
	0.50 mm	120 cm	RA66052M12	+CC	+22	480.00	+21	640.00	+23	580.00
	0.50 mm	130 cm	RA66052M13	+CC	+22	500.00	+21	660.00	+23	600.00
Size XXL	0.25 mm	100 cm	RA66022L10	+CC	+22	340.00	+21	430.00	+23	400.00
	0.25 mm	110 cm	RA66022L11	+CC	+22	360.00	+21	450.00	+23	420.00
	0.25 mm	120 cm	RA66022L12	+CC	+22	380.00	+21	470.00	+23	440.00
	0.25 mm	130 cm	RA66022L13	+CC	+22	400.00	+21	490.00	+23	460.00
	0.35 mm	100 cm	RA66032L10	+CC	+22	400.00	+28	470.00	+23	450.00
	0.35 mm	110 cm	RA66032L11	+CC	+22	420.00	+28	490.00	+23	470.00
	0.35 mm	120 cm	RA66032L12	+CC	+22	440.00	+28	510.00	+23	490.00
	0.35 mm	130 cm	RA66032L13	+CC	+22	460.00	+28	530.00	+23	510.00
Size XXXL	0.50 mm	100 cm	RA66052L10	+CC	+22	450.00	+21	610.00	+23	550.00
	0.50 mm	110 cm	RA66052L11	+CC	+22	470.00	+21	630.00	+23	570.00
	0.50 mm	120 cm	RA66052L12	+CC	+22	490.00	+21	650.00	+23	590.00
	0.50 mm	130 cm	RA66052L13	+CC	+22	510.00	+21	670.00	+23	610.00
	0.25 mm	100 cm	RA66022X10	+CC	+22	350.00	+21	440.00	+23	410.00
	0.25 mm	110 cm	RA66022X11	+CC	+22	370.00	+21	460.00	+23	430.00
	0.25 mm	120 cm	RA66022X12	+CC	+22	390.00	+21	480.00	+23	450.00
	0.25 mm	130 cm	RA66022X13	+CC	+22	410.00	+21	500.00	+23	470.00
Size XXXXL	0.35 mm	100 cm	RA66032X10	+CC	+22	410.00	+28	480.00	+23	460.00
	0.35 mm	110 cm	RA66032X11	+CC	+22	430.00	+28	500.00	+23	480.00
	0.35 mm	120 cm	RA66032X12	+CC	+22	450.00	+28	520.00	+23	500.00
	0.35 mm	130 cm	RA66032X13	+CC	+22	470.00	+28	540.00	+23	520.00
	0.50 mm	100 cm	RA66052X10	+CC	+22	460.00	+21	620.00	+23	560.00
	0.50 mm	110 cm	RA66052X11	+CC	+22	480.00	+21	640.00	+23	580.00
	0.50 mm	120 cm	RA66052X12	+CC	+22	500.00	+21	660.00	+23	600.00
	0.50 mm	130 cm	RA66052X13	+CC	+22	520.00	+21	680.00	+23	620.00

HS

DR

RP

RD

AW

AZ