

WURZELHEBER

Root elevators
 Elévateurs de racines
 Elevadores para raíces
 Leve per radice



KAPITEL

Chapter
 Chapitre
 Capítulo
 Capitolo

9



LINDO LEVIEN



4801 2,2 mm
4932 3,0 mm
6052 4,0 mm

Wurzelheber LINDO LEVIEN, gezahnt
Root elevator LINDO LEVIEN, with serration
Elévateur de racines LINDO LEVIEN, dentelé
Elevador de raíces LINDO LEVIEN, dentado
Leva per radice secondo LINDO LEVIEN, dentata

7226 4 mm links left gauche izquierdo sinistro




7227 4 mm rechts right droit derecho destro



Zahnung hat Widerhaken-Funktion – keine Hebelbewegungen!
 Serration for a barb effect – no leverage movement!

Les dents fonctionnent par petits mouvements latéraux comme un barbillon – ne pas faire un mouvement de levier!
 ¡Los dientes agarran por pequeños movimientos laterales, como un barbillon – no hacer movimientos de levantar!
 I denti ad uncino funzionano con piccoli movimenti laterali – non fare alcun movimento per sollevare!



HYLIN
 3,5 mm
 4243



HYLIN
 4,5 mm
 4244



HYLIN 30
 3,5 mm
 5646



HYLIN L
 3,5 mm
 4521



HYLIN R
 3,5 mm
 4522