



USC 128 Scaler	
#6	SUSC1286
#24	SUSC128



#1 Springfield Scaler, Single-End	
#6	SSP16
#24	SSP1



Younger-Good 15 Scaler, Single-End	
#6	SYG156
#1	SYG15



#WH-2 Whiteside Scaler, Single-End	
#24	SWH2

Posterior Sickle Scalers


Posterior Sickle Scalers are designed with contra-angled shanks for areas difficult to access, such as posterior interproximal surfaces. The blades can be curved or straight.




Learn more about our [Sickle Scalers](#).



204 Posterior Sickle Scaler	
#6	S2046
#2	S204



204S Posterior Sickle Scaler	
Harmony	S204SXE2
#9	S204S9E2
#G9, Steel Grip	S204SG9E2
#C8, Red	S204SC8E2
#6	S204S6
#7	S204S7
#8, Gray	S204S8
#2	S204S



204SD Posterior Sickle Scaler	
Harmony	S204SDXE2
#9	S204SD9E2
#G9, Steel Grip	S204SDG9E2
#C8, Red	S204SDC8E2
#6	S204SD6
#7	S204SD7
#8, Gray	S204SD8
#2	S204SD2