

DISPOSABLE VENTED (SLOTTED) ORAL EVACUATORS



- Non toxic phthalates free plastic
- Vented to reduce unnecessary tissue grab
- Soft rigid tubing, allows for greater patient comfort yet is strong enough for cheek retraction
- Thoroughly polished surface for smooth and quiet operation
- Fits standard 0.43" (11 mm) suction holders
- Colours: White or assorted (assorted only available in 1000's pack)

Standard Pack, 100 pcs/pack
 001 Disposable White **£3.15**
Bulk Pack, 1000 pcs/box
 001-B Disposable White **£29.00**
 001A-B Disposable Assorted Colours .. **£31.50**

DISPOSABLE VENTED/NON-VENTED 'S' CUT ORAL EVACUATORS



- Non toxic phthalates free plastic
- Vented and non-vented combination
- Soft rigid tubing, allows for greater patient comfort yet is strong enough for cheek retraction
- Fits standard 0.43" (11 mm) suction holders
- Colours: White or assorted (assorted only available in 1000's pack)

Standard Pack, 100 pcs/pack
 301 Disposable White **£3.15**
Bulk Pack, 1000 pcs/box
 301-B Disposable White **£28.50**
 301A-B Disposable Assorted Colours .. **£30.50**



Standard Pack, 100 pcs/pack
 300W White Tips/Clear Tubes **£2.75**
 300B Blue Tips/Clear Tubes **£2.75**
Saliva Ejectors Bulk Pack, 10 x 100 pcs/pack
 300W-B White Tips/Clear Tubes **£25.00**
 300B-B Blue Tips/Clear Tubes **£25.00**
 300A-B White Tips/Coloured Tubes ... **£27.00**

VENTED 'S' CUT ORAL EVACUATORS



- Non toxic phthalates free plastic
- Vented to reduce unnecessary tissue grab
- Soft rigid tubing, allows for greater patient comfort yet is strong enough for cheek retraction
- Thoroughly polished surface for smooth and quiet operation
- Fits standard 0.43" (11 mm) suction holders
- Colours: White or assorted (assorted only available in 1000's pack)


Standard Pack, 100 pcs/pack
 308 Disposable White **£3.15**
Bulk Pack, 1000 pcs/box
 308-B Disposable White **£29.00**
 308A-B Disposable Assorted Colours .. **£31.50**



DISPOSABLE SOFT TOUCH™ ORAL EVACUATORS (SLOTTED) VENTED LATEX FREE

- Non toxic phthalates free plastic
- Single use
- Soft Touch edge eliminates the risk of lacerating the oral tissue
- Rigid tubing strong enough for cheek retraction
- Vented to reduce unnecessary tissue grab and increase evacuation cover area
- Thoroughly polished surface for smooth and quiet operation
- Fits standard 0.43" (11 mm) suction holders
- Colours: White (Blue or Lavender Tip)

ST001 WHITE
 Soft Touch™ Oral Evacuators
 50 pcs/pack **£4.45**



AUTOCLAVABLE SURGICAL ASPIRATOR TIPS

- Non toxic phthalates free plastic
- Non-sterile, single use
- Autoclavable up to 134°C
- Fits standard 0.43" (11 mm) suction holders
- Colours: Standard - White Large orifice - Blue

Autoclavable Surgical Tips, 25 pcs/pack
 007 Mini Tip **£3.65**
 009 Standard Non-Vented **£4.15**
 009V Standard Vented **£5.15**
 008 Large Orifice Non-Vented .. **£4.15**
 008V Large Orifice Vented **£5.15**

AUTOCLAVABLE HVE SUCTION TUBES

- Fits European style 16 mm or standard 11 mm suction tubing
- Single use
- Autoclavable up to 134°C
- Special high quality plastic, sturdy enough for cheek retraction
- Fully polished internal surface and edge for quiet operation and patient comfort
- Assortment of bright, lively colours: Lavender, Light Blue, Light Green, Lime and Yellow



Assorted Colours, 10 pcs/pack
 335L Adult Large (16 mm) **£7.45**
 335M Adult Medium (16 mm) .. **£7.45**
 336P Pedo (16 mm) **£7.45**

Assorted Colours, 10 pcs/pack
 005L Adult Large (16 mm) **£7.45**
 005M Adult Medium (16 mm) .. **£7.45**
 006P Pedo (16 mm) **£7.45**

Assorted Colours, 10 pcs/pack
 205M Adult Medium (11 mm) .. **£7.45**
 206P Pedo (11 mm) **£7.45**

AUTOCLAVABLE SUCTION TIP SLEEVES




- For protection from sharp edges on suction tubes
- Autoclavable up to 134°C
- Colour: Light Blue

290 Suction Tip Sleeves 6 pcs/pack **£9.95**

SUCTION TUBE ADAPTORS

- Autoclavable up to 134°C

Buy 2 get 1 FREE
 (same item)




Saliva Ejector 7 mm
 1100 Saliva Ejector Adaptor reduces 11 mm to 7 mm
 Standard HVE Hose 11 mm **£12.95**

Saliva Ejector 6.5 mm
 306 Saliva Ejector Adaptor reduces 11 mm to 6.5 mm
 Standard HVE Hose 11 mm **£12.95**

Suction Tube 11 mm
 305 Suction Tube Adaptor reduces 16 mm to 11 mm
 European HVE Hose 16 mm **£2.95**

ASPIRATOR BRUSHES

- Autoclavable up to 134°C
- Rigid bristles for effective cleaning



6 pcs/pack
 020 X-Fine dia. 5 mm **£7.45**
 021 Fine dia. 8 mm **£7.95**
 022 Medium dia. 14 mm .. **£7.95**
 023 Large dia. 20 mm **£8.95**