

Endo-Bite™

X-ray holders for diagnostic X-rays

The Endo-Bite™ is a special film and phosphor plate holder for diagnostic X-rays of root canal instruments during endodontic treatment, checks on root fillings, as well as inserted pins and screws, and for general use when the bite is blocked (e.g. rubber dam clamps). The system includes two types of holder – anterior (green) and posterior (red) – for use in all four quadrants.

Compatibility

The Endo-Bite is also suitable for indirect digital systems. It is compatible with all major brands of phosphor plate sensors.

Endo-Bite with Ring Assorted Kit

REF/Part No. 1790

2 Endo-Bite Anterior with Ring, 2 Endo-Bite Posterior with Ring, 2 Centring Devices

Endo-Bite with Ring Refill Packs

REF/Part No.

4 Endo-Bite Anterior with Ring, 1 Centring Device


1791

4 Endo-Bite Posterior with Ring, 1 Centring Device

1792



Endo-Bite Anterior



Endo-Bite Posterior

Paro-Bite™

X-ray for perio patients and children

The Paro-Bite™ is a special film holder for assessing the bone level in perio patients and support zones in children (vertical bite wing). It guarantees radiographs without overlaps and therefore eliminates the need for repetition and unnecessary exposure to X-rays. A direct comparison of radiographs over the years can be made for the same patient.


Compatibility

The Paro-Bite™ is compatible with all major brands of phosphor plate sensors.

Paro-Bite with Ring

REF/Part No. 1770

5 Paro-Bite with Ring, 1 Centring Device



Paro-Bite with Ring