





RUBBER DAM CLAMPS


UPPER MOLARS




4
| RDCM4
Flat narrow jaws for small upper molars.



W4
| RDCMW4
Flat narrow jaws for small upper molars.



5
| RDCM5
Slightly festooned jaws for large upper molars.




8
| RDCM8
Deeply festooned jaws for upper molars.




26N
| RDCM26N
Flat jaws for upper molars.




30
| RDCM30
Offset jaws for upper right molars; paired with RDCM31.



31
| RDCM31
Offset jaws for upper left molars; paired with RDCM30.




201
| RDCM201
Slightly festooned jaws for upper molars.




205
| RDCM205
Slightly festooned jaws for large upper molars.


LOWER MOLARS



3
| RDCM3
Flat narrow jaws for small lower molars.



W3
| RDCMW3
Flat narrow jaws for small lower molars.




7
| RDCM7
Flat jaws for lower molars.




7A
| RDCM7A
Slightly reversed festooned jaws for large lower molars.



W7
| RDCMW7
Flat jaws for lower molars.



12A
| RDCM12A
Serrated offset jaws for lower right molars; paired with RDCM13A.



13A
| RDCM13A
Serrated offset jaws for lower left molars; paired with RDCM12A.



26
| RDCM26
Slightly reversed festooned jaws for small lower molars.



28
| RDCM28
Flat jaws for lower molars.