

Bite Registration Material



Alminax Bite Wax

Key Features

Colour: Silver Grey
 Sheet size: 143mm x 73mm x 2.6mm
 Softens at 55°C
 A proven bite registration material
 Dense yet easily softened bite wax, stable at room temperature

Packaging: 250g/10 sheets, 500g/20 sheets



Alminax Wax Full Arch

Key Features

Colour: Silver Grey
 Thickness: 2.6mm
 Softens at 55°C
 Ready shaped to mouth, proven bite registration material, dense yet easily softened bite wax, stable at room temperature

Packaging: 20 per box, 40 per box



Rite Bite Index

Key Features

Colour: Silver Grey
 Thickness: 2.6mm
 Softens at 55°C
 Ready shaped to mouth, reinforced around inner edge. A proven bite registration material, dense yet easily softened bite wax, stable at room temperature

Packaging: 30 per box/24 reinforced + 6 plain



Ortho Bite Wax

Key Features

Colour: Deep Red
 Sheet size: 142mm x 73mm x 2.2mm
 Soft pliable wax specially developed for bite registration in orthodontic cases

Packaging: 500g box



Detex Articulating Paper

Key Features

Colour: Blue
 Sheet size: 120mm x 25mm x 40 micron thick, double sided
 Special vegetable wax coating, pure colouring ingredient does not stain lips
 Provides clear sharp mark on natural or artificial teeth to show where bite adjustment is necessary
 Stitched interleaving between sheets keeps paper in perfect condition until needed and prevents curling.

Packaging: 12 books per box

Cross Infection Control Products



PracticeSafe Disinfectant

Key Features

Disinfectant for hard alcohol resistant surfaces in Dental Practices and Commercial Laboratories plus medical products such as hand and angle pieces, impressions, try-ins and models
 Free from: aldehydes, phenols, sulfates, azoles and peroxides
 Effective Spectrum:
Bacteria: TB, E. Coli, Enterococcus hirae, Staphylococcus aureus, Pseudomonas aeruginosa
Viruses: Vaccinia, Hepatitis B, Hepatitis C, HIV, BVDV
Fungus: Fungicidal
 Fragrance: Lemon, Fruit, Flower, Natural
 No overpowering odour

Packaging: 500ml, 1L, 5L refill, 10L refill

Fragrances



PracticeSafe Disinfectant Wipes

Key Features

Available in 100 Microfibre Wipes, 250 Heavy Duty Microfibre Wipes or 600 Economy Wipes complete tubs.
 100 Microfibre Wipes are refillable using the 100 Microfibre complete refill.
 Prepared using PracticeSafe Disinfectant for hard alcohol resistant surfaces in Dental Practices, Commercial Laboratories, medical products such as hand and angle pieces, impressions, try-ins and models
 Free from: aldehydes, phenols, sulfates, azoles and peroxides

Effective Spectrum:
Bacteria: TB, E. Coli, Enterococcus hirae, Staphylococcus aureus, Pseudomonas aeruginosa
Viruses: Vaccinia, Hepatitis B, Hepatitis C, HIV, BVDV
Fungus: Fungicidal
 Sheet size: 100 Microfibre 18.5cm x 13.5cm, 250 HD Microfibre 26cm x 17cm, 600 Economy 24cm x 14.5cm
 Fragrance: Lemon, Fruit, Flower, Natural
 No overpowering odour, non-drip, durable and gentle on hands

Packaging: 100 Microfibre sheets, 100 Microfibre Sheet Refill, 250 HD Microfibre sheets, 600 Economy sheets

Fragrances



Heavy Duty

Economy

Heavy Duty Microfibre & Economy Dry Rolls

Key Features

The Dry Rolls can be used with ChairSafe Alcohol Free Disinfectant or PracticeSafe Disinfectant to disinfect surfaces

Packaging: HD Refill 2 x 250 sheets, Economy Refill 2 x 600 Sheets