

Universal Tray Adhesive

Adhesive for A-Silicones

Impression systems / Accessories

Universal Tray Adhesive is a specific adhesive for A-Silicones.

Universal Tray Adhesive is an adhesive designed to ensure the adhesion of addition silicone to standard metal solid impression trays and individual resin plates.

Characteristics

- Brushing adhesive
- Blue colour

Advantages

- Prevents the separation between addition silicones and metal or resin trays, by reducing the risk of distortion during the removal of the impression from the oral cavity.
- Fast solvent evaporation
- Easily removable residues



Universal Tray Adhesive

Code	Packaging
C700025	10 ml bottle