

VOCO GIC Mixer®

Capsule mixing device especially for glass ionomer cements



REF 9050 Capsule mixing device

Indications

Device for mixing powder and liquid in application capsules

Advantages

- High oscillating frequency of 4,600 rpm, creating a particularly homogeneous mixture of powder and liquid
- Modern design
- Capsule holder allows simple insertion and removal of capsules
- User-friendly
- Quiet and practically vibration-free
- Mixing times can be freely selected between 1 and 99 seconds

VOCO Dynamic Dispenser

Dynamic mixing device



REF 9140 Mixing device 230 V
Not available in Canada

Indications

Device for the automatic dispensing and mixing of VOCO impression materials in the 5:1 cartridge

Advantages

- Homogeneous, dynamic mixing of product components to ensure a consistently high material quality
- Bubble-free dispensing
- Fill level display
- 2 feed speeds
- Electronic residue detection