



#TU-17/23 Wilkins Tufts DE Explorer

#6	EXTU17/236
#7	EXTU17/237
#8	EXTU17/238
#31	EXTU17/23



Combination #23 Shepherd's Hook and #6 DE Explorer

#6	EX23/66
#31	EX23/6



#5 DE Explorer

#6	EXD56
#7	EXD57
#8	EXD58
#31	EXD5



Combination #6 and #3 DE Explorer

#6	EXD66
#31	EXD6



Explorer 3A/TU17

#6	EX3A/TU176
#31	EX3A/TU17




#11/12 ODU DE Explorer

#6	EXD11/126
#7	EXD11/127
#8	EXD11/128
#31	EXD11/12



#11/12 ODU After Five DE Explorer

#6	EXD11/12A6
#7	EXD11/12A7
#8	EXD11/12A8
#31	EXD11/12AF



To correctly adapt the tip of the 11/12 ODU explorer, place the tip interproximally so that the terminal shank is parallel to the tooth surface.

XTS™ Explorers

XTS Explorers feature a non-stick Aluminum Titanium Nitride coating. It is perfect for contouring dental anatomy during composite restorations. The dark coating also improves contrast between tooth and instrument tip.



#2 DE XTS™ Explorer

#6	EXD2XTS
----	---------




#3 Cowhorn DE XTS™ Explorer

#6	EXD3CHXTS
----	-----------



#5 DE XTS™ Explorer

#6	EXD5XTS
----	---------



#TU17/23 DE XTS™ Explorer

#6	EXTU1723XTS
----	-------------