





Jacquette Sickle Scalers

Jacquette Sickle Scalers have two cutting edges on a straight blade that end at a sharp point. The back surface of the blade is pointed giving a triangular cross-section.




The Curved Sickle Scaler
Two cutting edges on a curved blade that end in a sharp point




The Straight Sickle Scaler
Two cutting edges on a straight blade that end in a sharp point (also known as a Jacquette Scaler)



30/33 Standard/Mini Jacquette Scaler	
Harmony	SJ30/33XE2
#9	SJ30/339E2
#C8, Lavender	SJ30/33C8E2
#6	SJ30/336
#7	SJ30/337
#8, Gray	SJ30/338
#2	SJ30/33



31/32 Standard Jacquette Scaler	
Harmony	SJ31/32XE2
#9	SJ31/329E2
#C8, Red	SJ31/32C8E2
#6	SJ31/326
#7	SJ31/327
#8, Gray	SJ31/328
#2	SJ31/32



34/35 Mini Jacquette Scaler	
Harmony	SJ34/35XE2
#9	SJ34/359E2
#C8, Red	SJ34/35C8E2
#6	SJ34/356
#7	SJ34/357
#8, Gray	SJ34/358
#2	SJ34/35



H3/H4 Jacquette Scaler	
#6	SH3/46
#2	SH3/4



14/15 Posterior Jacquette Scaler	
#6	SJ14/156
#2	SJ14/15