



Charisma® Flow
Naturally reliable.

Giving a hand to oral health.



KULZER
MITSUI CHEMICALS GROUP

Charisma® Flow – The flowable all-rounder.

Charisma Flow is the ideal supplement to the line of universal Charisma composites. Its low viscosity and good thixotropy makes it the ideal choice for many indications, such as cavity lining, fissure sealing and minimal invasive fillings.

Charisma Flow stands for safe usage and a high standard of naturally reliable restorations

- Versatile all-rounder for many indications.
- Easy and fast handling due to its low viscosity and good thixotropy.
- Its radiopacity enables a safe and reliable diagnosis.
- Pronounced chameleon effect.
- Shades of Charisma Flow optimally match the Charisma shades.

The **Charisma Flow Baseline** is a white and opaque shade that can be applied on the cavity floor to mark the transition to the natural hard tooth tissue. This means that unnecessary drilling away of healthy hard tooth tissue when removing an insufficient filling can be avoided.



The white opaque Charisma Flow Baseline marks the cavity floor.
Photograph: Prof. Dr. C.-P. Ernst, University of Mainz, Germany

Art. code	Article
66015521	Charisma Flow Syringe Assortment 1 x 1.8 g Syringe each of the shades: A1, A2, A3, Baseline



Shades	Charisma Flow
	Syringe 1 x 1.8 g
Universal shades	
A1	66015495
A2	66015496
A3	66015497
A3.5	66015498
B2	66015499
B3	66015500
Opaque shades	
OA2	66015501
OA4	66015502
Special shade	
Baseline	66015504