

DIRECTA PrepSure

PrepSure is a set of three innovative crown preparation instruments designed to help you monitor and confirm the amount of tooth reduction for perfect occlusal clearance. With three dedicated gauges, PrepSure takes your crown preparations to the next level, ensuring greater accuracy, confidence and efficiency every time.


These instruments are essential for improving the accuracy and efficiency of crown preparation.

- Takes out all the guesswork.
- Avoids under or over-reduction of tooth structure.
- Prevents a second impression or digital scan.
- Solves a common problem with crown preparation.
- A suitable complement to your digital workflow.
- Autoclavable.

Achieve fast and reliable tooth preparation with PrepSure, an instrument that ensures accurate occlusal clearance for crowns and other prosthetic restorations.

PrepSure ensures that every preparation meets the specific requirements of the material used, minimising errors and enhancing clinical outcomes. PrepSure is durable, precise and essential for modern crown workflows, helping you to avoid under- or over-reduction across all restorative materials.

PrepSure, Crown Prep Guide



M. Tip confirms mesial-half clearance




D. Tip confirms distal-half clearance


DIRECTA PrepSure

KEY FEATURES


- Indicates and confirms the space for tooth reduction required for crowns and other prosthetic restorations.
- Compatible with multiple prosthetic materials.




Prod No **Product Description**
92410 PrepSure Single 1.0 mm green, 1 pc
92413 PrepSure 1.0 mm green Set, 3 pcs




Prod No **Product Description**
92415 PrepSure Single 1.5 mm red, 1 pc
92418 PrepSure 1.5mm red Set, 3 pcs



Prod No **Product Description**
92420 PrepSure Single 2.0 mm blue, 1 pc
92423 PrepSure 2.0 mm blue Set, 3 pcs



Prod No **Product Description**
92403 PrepSure Assorted Set
(1.0 mm, 1.5 mm, 2.0 mm), 3 pcs



For more information about PrepSure
scan the QR-code with your mobile
or visit www.directadental.com