

POLO|DENT[®]

DENTAL DAM

Dental Dam

152x152 MM



Product information

Issued 12/2023

Dental Dam

Made from natural rubber latex or latex free

PoloDent dental dams are used within dentistry to isolate specific areas of the teeth and gums during procedures. They improve infection control, visibility and cleanliness of the working area and patient experience. Their use is now acknowledged around the world as best practice for all types of dentistry.

PoloDent Dental Dam is available in two main material options based on natural rubber latex and latex free (synthetic natural rubber). The formulations for both materials are optimised to give high tear resistance and the right elastomeric properties. They are also produced using FDA approved ingredients. Low protein latex is used to minimise any risk of sensitivity. Our latex free dams use **a specially synthesised version of the natural rubber molecule to closely match the physical properties of natural rubber with no proteins.**



Dental Dam Blue Medium NON-LATEX

thickness 0,18-0,23 mm
contains 20 pieces
REF 34 39 00

Dental Dam Blue Heavy NON-LATEX

thickness 0,23-0,28 mm
contains 20 pieces
REF 34 39 05

Dental Dam Green Medium LATEX

thickness 0,18-0,23 mm
contains 36 pieces
REF 34 39 10

Dental Dam Green Heavy LATEX

thickness 0,23-0,28 mm
contains 36 pieces
REF 34 39 15