



Mouth Props



Mouth Prop, Latex-Free

Pair of 2

Small Child [MP48](#)

Child [MP52](#)

Adult [MP54](#)

Needle Holders

Our Perma Sharp™ Needle Holders feature tungsten carbide inserts with a matched diamond-cut pyramidal profile for enhanced longevity, a secure, positive suture needle grasp and reduced risk of cutting the suture material.



Learn more about our [Needle Holders](#).

Microsurgical Needle Holders

Microsurgical needle holders have a lock to ensure that the needle is firmly secured. This allows the needle to be guided through coarse gingival tissue with controlled grip pressure.



Castroviejo Microsurgical Needle Holder, Straight

5-0, 6-0, 7-0, 8-0

14 cm

[NH5020M](#)



Castroviejo Microsurgical Needle Holder, Curved

5-0, 6-0, 7-0, 8-0

14 cm

[NH5021M](#)