



Probes



Qulix Probes are coded with accurate, smooth black markings, which will not chip, flake, or fade. Colorvue™ Probes feature a vivid markings for enhanced contrast intraorally.



Learn more about our [Probes](#).

Colorvue™ Probes

Colorvue probes feature a vivid yellow or white tip with black markings for enhanced contrast to intraoral structures in both natural and implant dentitions. Optimal flexibility and a rounded tip provide greater patient comfort and acceptability. The probe tips are made with resin, which makes them safe for use with implants and feature a twist-on design for easy replacement of worn tips.



Colorvue is a safe and effective choice for assessing peri-implant health.²



#12 Colorvue™ Probe Tips

12 tips (3-6-9-12)

[PCV12PT](#)

#12 Colorvue™ Probe Kit

1 Handle + 7 tips
(3-6-9-12)

[PCV12KIT6](#)

2 Handles + 12 tips
(3-6-9-12)

[PCV12KIT12](#)



#11 Colorvue™ Probe Tips

12 tips (3-6-8-11)

[PCV11PT](#)

#11 Colorvue™ Probe Kit

1 Handle + 7 tips
(3-6-8-11)

[PCV11KIT6](#)

2 Handles + 12 tips
(3-6-8-11)

[PCV11KIT12](#)



#UNC12 Colorvue™ Probe Tips

12 tips

(1-2-3-4-5-6-7-8-9-10-11-12)

[PCVUNC12PT](#)

#UNC12 Colorvue™ Probe Kit

1 Handle + 7 tips

(1-2-3-4-5-6-7-8-9-10-11-12)

[PCVNCKIT6](#)

2 Handles + 12 tips

(1-2-3-4-5-6-7-8-9-10-11-12)

[PCVNCKIT12](#)



Colorvue™ Oxford Probe Tips

12 tips

(1-2-3-4-5-6-7-8-9-10-11-12-13-14-15)

[PCVO12PT](#)