



Upper Molar Clamps



Rubber Dam Clamp #4

Flat narrow jaws for small upper molars.

Black Line **RDCM4X**

Satin Steel™ Alloy **RDCM4**



Rubber Dam Clamp #5

Slightly festooned jaws for large upper molars.

Black Line **RDCM5X**

Satin Steel™ Alloy **RDCM5**



Rubber Dam Clamp #26N

Flat jaws for upper molars.

Black Line **RDCM26NX**

Satin Steel™ Alloy **RDCM26N**



Rubber Dam Clamp #31

Offset jaws for upper left molars; paired with RDCM30.

Black Line **RDCM31X**

Satin Steel™ Alloy **RDCM31**



Rubber Dam Clamp #205

Slightly festooned jaws for large upper molars.

Black Line **RDCM205X**

Satin Steel™ Alloy **RDCM205**



Rubber Dam Clamp #W4

Flat narrow jaws for small upper molars.

Black Line **RDCMW4X**

Satin Steel™ Alloy **RDCMW4**



Rubber Dam Clamp #8

Deeply festooned jaws for upper molars.

Black Line **RDCM8X**

Satin Steel™ Alloy **RDCM8**



Rubber Dam Clamp #30

Offset jaws for upper right molars; paired with RDCM31.

Black Line **RDCM30X**

Satin Steel™ Alloy **RDCM30**



Rubber Dam Clamp #201

Slightly festooned jaws for upper molars.

Black Line **RDCM201X**

Satin Steel™ Alloy **RDCM201**