



*Cavex laboratory materials are the natural extension of the dental practice and must therefore meet stringent quality requirements. Thanks to Cavex's close collaboration with dental technicians and prosthetists these materials are also closely matched to the specific requirements of the dental laboratory. For dental use only.*

## **How nature and Cavex provide a stable basis for perfect prostheses.**

Shellac is one of the natural ingredients of Cavex Dental Baseplates and Cavex Individual Trays. The modern Cavex technology makes it possible to reshape a two-dimensional plate into an extremely stable 3D object. Simply through the application of heat and cooling down.

For you it's the most reliable basis for perfectly-fitting prostheses.



## Cavex Dental Baseplates perfected: Mica-free.

- **rapid processing time**
- **dimensionally stable**
- **economical**
- **standard thickness**

The pink plate, standard thickness 1.4 mm, can easily be shaped to the gypsum model after brief heat penetration in the flame, and is easy to cut. After curing it will assume its final form, which will not spring back. Mouth temperature no longer has any effect. A stable basis for the wax/bite wafer.

The latest product development - Mica-free - means a further improvement in the performance of Cavex Dental Baseplates. The presence of aluminium hydroxide ensures that - above 200°C - cooling water of crystallisation is released, which prevents (or at least retards) burning of the plate. The end result is also smoother, non-porous and has an attractive sheen, just like the later prosthetic plastic.

## Cavex Individual Trays: natural reproduction that remains rigid.

The proven effectiveness of the grey-coloured impression tray, standard thickness 1.9 mm, guarantees an optimum fit, even if the tray is stored for quite a long period. Its dimensional stability is proverbial, and the combination of a rapid processing time and a favourable price/performance ratio has ensured that Cavex Individual Trays have already been popular in dental laboratories for many years.

The Dental Baseplates, Uppers and Lower, and the Individual Trays, Uppers and Lower, are packaged in handy, robust boxes of 100 pieces. The metal handles, which are easy to fix into the still warm material, are included as standard.

Cavex Dental Baseplates and Cavex Individual Trays are available from your supplier. For further information on product composition and applications, please contact Cavex, +31 23 530 77 00 - Dental Customer Service, or surf to [www.cavex.nl](http://www.cavex.nl)

