

Periotomes

Our Periotomes are designed to facilitate the removal of teeth with minimal damage to the surrounding alveolar bone by severing the periodontal ligament with the use of the sharp, thin clinically designed blades.



Learn more about our [Periotomes](#).



Periotome, Posterior

Black Line	PT1X
#6	PT1



Periotome, Anterior

Black Line	PT2X
#6	PT2

Periotome, Posterior, Right

Black Line	PT3X
#6	PT3



Periotome, Anterior, Angled

Black Line	PT5X
#6	PT5

Periotome, Posterior, Left

Black Line	PT4X
#6	PT4



Periotome, Anterior, Straight

Black Line	PT6X
#6	PT6



Wise Periotome

Black Line	PEREURX
#6	PEREUR