



Universal Currettes

Universal Currettes have two parallel cutting edges that meet at a rounded toe. The back surface is rounded giving a semi-circular cross section. Due to this rounded back surface, currettes are ideal for the removal of subgingival debris.



Learn more about our [Universal Currettes](#).

Langer Currettes

The Langer Curette design combines the benefits of a universal curette blade with two cutting edges and the shank designs of area specific currettes. The Langer Currettes were originally designed to be used in a set of three (1/2, 3/4, 5/6), to adapt to most tooth surfaces, and to be suitable for scaling as well as periodontal surgery.



Langer 1/2 Curette	Standard	Rigid
Harmony	SL1/2XE2	N/A
#9	SL1/29E2	N/A
#C8, Lavender	SL1/2C8E2	N/A
#6	SL1/26	SL1/2R6
#7	SL1/27	SL1/2R7
#8, Gray	SL1/28	N/A
#2	SL1/22	N/A

Langer 1/2 After Five™ Curette	
A 3 mm longer terminal shank and a thinner blade for improved access	
#9	SL1/2AF9E2
#6	SL1/2AF6

Langer 1/2 Mini Five™ Curette	
Longer terminal shank, standard diameter, and blade half the length of the After Five™ Currettes	
#9	SL1/2MF9E2
#6	SL1/2MF6

Best in Practice Fact

- The Langer curette takes a “best of both worlds” approach to instrument design. The curette pairs a universal blade — featuring two parallel edges and a 90-degree angle between the terminal shank and blade face — with the shank design of area-specific Gracey currettes. The result is an instrument that adapts well and is efficient.