



# TOP TRAY

## DISPOSABLE TRAY PAPER

### CONFIGURATION OVERVIEW

**Type:**

- Non-Sterile
- Wood free
- Latex free

**Material:**

- Natural cellulose fibers

### QUALITY SPECIFICATION

**Conformity to standards:**

- DIN EN ISO 536,
- DIN EN 20 534,
- ISO 534,
- DIN EN ISO 1924-2,
- DIM ISO 3781

### PRODUCT VARIATIONS

<b>Color code</b>	WR	L	Y	P	F	F
<b>Color</b>	Burgundy	Lilac	Yellow	Pink	Fuchsia	Green
						
<b>Color code</b>	FG	B	MB	W	O	
<b>Color</b>	Fresh Green	Blue	Med Blue	White	Orange	
						

## TECHNICAL SPECIFICATIONS

### Optical brightener:

- negative

### Chlorine:

- negative

Parameter	Standard	Data
Basis weight	DIM EN ISO 536	55 ± 2 g/m <sup>2</sup>
Thickness	DIM EN 20 534, ISO 534	0.15 ± 0.02 mm
Tensile strength MD	DIM EN ISO 1924-2	≥ 28 N/15 mm
Tensile strength CD	DIM EN ISO 1924-2	≥ 12 N/15 mm
Tensile strength MD, wet	DIM ISO 3781	≥ 4.0 N/15 mm
Tensile strength CD, wet	DIM ISO 3781	≥ 2.0 N/15 mm
Stretch MD	DIM EN ISO 1924-2	≥ 4 %

## PACKAGING

<b>Sheet/unit:</b>	250	<b>Units/carton:</b>	10 (size 18x28 cm), 5 (size (28x36 cm)
<b>Carton size:</b>	37.5x29.5x23 cm	<b>Carton weight:</b>	7.3 kg
<b>Carton/layer:</b>	8	<b>Layers/pallet:</b>	8
<b>Pallets/carton:</b>	64	<b>Pallet weight:</b>	487 kg
<b>Pallet size:</b>	80x120x2198 cm		