



endosuccess retreatment



ET20
F88011



ET20D
F88013



ET25
F88018



ET25S
F88021

Canal Retreatment



Retreatment tip, length 20mm, taper 6%

Used in the 1st coronal third:

- Extraction of filling material, silver points, broken instruments.
- Removal of debris and the smear layer.



Diamond-coated retreatment tip, 30 µm, length 20mm, taper 5%

Used in the 1st coronal third to remove very hard materials by brushing the walls.

The diamond coating of the ET20D tip increases the cutting and lateral abrasion effect.



Titanium-Niobium tip, length 20mm, taper 3%

Retreatment in the middle and apical thirds and the extraction of broken instruments.

The Titanium-Niobium alloy of the ET25 range allows perfect transmission of the ultrasonic vibrations and tip flexibility*.



Short Titanium-Niobium tip, length 15mm, taper 4%

Retreatment in the coronal third and the isthmuses.