



Langer 3/4 Curette	Standard	Rigid
Harmony	SL3/4XE2	N/A
#9	SL3/49E2	N/A
#C8, Lavender	SL3/4C8E2	N/A
#6	SL3/46	SL3/4R6
#7	SL3/47	SL3/4R7
#8, Gray	SL3/48	N/A
#2	SL3/42	N/A

Langer 3/4 Mini Five™ Curette

Longer terminal shank, standard diameter, and blade half the length of the After Five™ Curettes

#9	SL3/4MF9E2
#6	SL3/4MF6

Langer 3/4 After Five™ Curette

#9	SL3/4AF9E2
----	----------------------------



Langer 5/6 Curette	Standard	Rigid
Harmony	SL5/6XE2	N/A
#9	SL5/69E2	N/A
#C8, Lavender	SL5/6C8E2	N/A
#6	SL5/66	SL5/6R6
#7	SL5/67	N/A
#8, Gray	SL5/68	N/A

Langer 5/6 Mini Five™ Curette

Longer terminal shank, standard diameter, and blade half the length of the After Five™ Curettes

#9	SL5/6MF9E2
#6	SL5/6MF6



Langer 5/6 After Five™ Curette

#6	SL5/6AF6
----	--------------------------



Langer 17/18 Curette	Standard	Rigid
Harmony	SL17/18XE2	N/A
#9	SL17/189E2	N/A
#C8, Lavender	SL17/18C8E2	N/A
#6	SL17/186	SL17/18R6
#7	SL17/187	N/A
#8, Gray	SL17/188	N/A

Langer Application Chart



- Langer After Five and Mini Five Curette terminal shanks are elongated by 3 mm for access into deep periodontal pockets and root surfaces of 5 mm or more.
- Blade is 10% thinner than standard Langer to ease gingival insertion and reduce tissue distention.
- Langer Mini Five Curette has a reduced blade, half the length of the After Five Langer Curette or standard Langer, for better adaptation in narrow pockets and furcations.