



Hartmetall-Finierer


Carbide finishers · Fraises à finir au carbure de tungstène · Fresas para acabar

C48


ID-No.	009	010	012	014	016
500 104 243 071 ...	009	010	012	014	016
500 204 243 071 ...	009	010	012	014	016
500 314 243 071 ...	009	010	012	014	016
L mm	3,7	4,0	4,3	4,7	4,9

C243 K

ID-No.	014	016
500 204 297 072 ...		016
500 314 297 072 ...	014	
L mm	6,0	6,0

C244


ID-No.	012	014
500 204 289 072 ...	012	
500 314 289 072 ...	012*	014*
314 *max. 120.000 min ⁻¹		
L mm	8,0	8,0

C244 K


ID-No.	016	021
500 204 298 072 ...		021
500 314 298 072 ...	016*	021*
314 *max. 120.000 min ⁻¹		
L mm	8,0	8,0

Kronentrenner

crown cutters · maxi cut métal · tagliacorone

C21 RX

ID-No.	010	012
500 314 137 019 ...	010*	012*
L mm	4,0	4,3
314 *max. 160.000 min ⁻¹		

C23 RX

ID-No.	010	012
500 314 194 019 ...	010*	012*
L mm	4,0	4,3
314 *max. 160.000 min ⁻¹		

**C21 RC**

ID-No.	012
500 314 138 008 ...	012*
L mm	3,0
314 *max. 160.000 min ⁻¹	

H · HP · PM · PM

... **104** ...

W · RA · CA · CA

... **204** ...

FG · FG · FG · FG

... **314** ...