

A large range of tips meeting all the clinical needs

periodontics
 Periodontal debridement, Root planing
 N° 1S, H3, H4L, H4R tips,
 4 autoclavable dynamometric wrenches



F00936

perio maintenance BDR
 Biofilm disruption
 TK1-1S, TK1-1L, TK2-1L, TK2-1R
 tips, 4 autoclavable dynamometric
 wrenches



F00737

scaling
 Supra- and sub-gingival scaling
 N° 1, N° 1S, N° 10X, H3 tips,
 4 autoclavable dynamometric
 wrenches



F00934

perio precision
 Periodontal maintenance
 P2L, P2R, TK1-1S tips, 3 autoclavable
 dynamometric wrenches



F00939

implantprotect pure titanium
 Treatment of peri-implantitis and
 maintenance
 IP1, IP2L, IP2R, IP3L, IP3R tips,
 autoclavable metal support and
 universal wrench



F02120

perfect margin rounded
 Prosthetic finishing with
 chamfered shape
 PM1, PM2, PM3, PM4 tips,
 autoclavable metal support and
 universal wrench



F00738

hygiene
 Versatile, gentle hygiene treatment
 N° 1, N° 1S, N° 10Z, TK1-1S tips,
 4 autoclavable dynamometric
 wrenches



F00935

endo one
 Endodontic treatments
 CAP1, CAP2, CAP3, ET25, ETPR
 tips, 4 Irrisafe 25-21 mm blister,
 autoclavable metal support
 and universal wrench



F00732

endosuccess retreatment
 Canal Retreatment
 ET18D, ET20, ET25, ET25S,
 ETBD, ETPR tips, autoclavable
 metal support and universal wrench



F00737

perfect margin shoulder
 Prosthetic finishing with
 shoulder shape
 PMS1, PMS2, PMS3, PM4 tips,
 autoclavable metal support and
 universal wrench



F00736

excavus
 Minimal excavation and micro-
 abrasion
 EX1, EX2, EX3, EX-L, EX-R tips,
 autoclavable metal support and
 universal wrench




F00739

endosuccess canal access prep
 Canal access preparation
 CAP1, CAP2, CAP3 tips,
 autoclavable metal support and
 universal wrench



F88180

endosuccess apical surgery
 Apical Surgery
 AS3D, AS6D, AS9D, ASLD, ASRD
 tips, autoclavable metal support and
 universal wrench



F00069

perfect margin veneers
 Ceramic veneers finishing
 PMV1, PMV2, PMV3, PMV4 tips,
 autoclavable metal support and
 universal wrench



F02020