

Fully equipped

Supplementary set for VistaScan

Our supplementary sets provide all the materials you need for using the VistaScan models.



Supplementary set VistaScan Combi View / Perio Plus

Image plate size 1, size 3, 2 pcs each

Image plate size 4, 1 pce

Foil cassette size 1, size 3, size 4, 1 pce each

Light protection cover size 1, size 3, size 4, 100 pcs each

Bite block protector size 4, 100 pcs

2130-007-00 VistaScan Combi View

2140-007-00 VistaScan Perio Plus

Supplementary set VistaScan Mini Plus /

VistaScan Mini View

Image plate size 1, size 3, 2 pcs each

Image plate size 4, 1 pce

Light protection cover size 1, size 3, size 4, 100 pcs each

Bite block protector size 4, 100 pcs

2141-007-00



Highest image resolution

Image plate plus

Thanks to an image resolution of 22 LP/mm, Image Plates Plus are superior to X-ray film when it comes to picking up detail – and therefore superior to the competition as well. The reusable image plates are particularly thin and flexible, and therefore comfortable for patients.



Our range of intraoral image plates

Image Plates Plus size 0

Format 2 x 3 cm, 2pcs

2130-040-50

2130-040-60 with copper dot and ID

Image Plates Plus size 1

Format 2 x 4 cm, 2 pcs

2130-041-50

Image Plates Plus size 2

Format 3 x 4 cm, 4 pcs

2130-042-50

2130-042-60 with copper dot and ID

Image Plates Plus size 3

Format 2.7 x 5.4 cm, 2 pcs

2130-043-50

Image Plates Plus size 4

Format 5.7 x 7.6 cm, 1 pcs

2130-044-50

PO07-597-02/IOD-dt.de/1/01/H3 Subject to technical changes

DÜRR DENTAL SE
Höfigheimer Str. 7
74321 Bietigheim-Bissingen
Germany
www.duerrendental.com
info@duerrdental.com

 **DÜRR
DENTAL**
THE BEST, BY DESIGN

VistaScan image plate accessories



COMPRESSED AIR
SUCTION
IMAGING
DENTAL CARE
HYGIENE

 **DÜRR
DENTAL**
THE BEST, BY DESIGN

Safe and reliable to use

Optimum protection

The light protection covers for Image Plates Plus for protection against incident light and contamination are made of material that has been certified to comply with health requirements, they are 100 % hygienically leak-tight and have narrow edges and a notch on both sides for left and right-handed users. Their exceptional impermeability to light has been demonstrated in studies.



Light protection cover size 0

Protective cover for image plate, format 2 x 3 cm, 100 pcs

2130-080-00/2130-080-50*

Light protection cover size 1

Protective cover for image plate, format 2 x 4 cm, 100 pcs

2130-081-00

Light protection cover size 2

Protective cover for image plate, format 3 x 4 cm, 300 pcs

2130-082-00/2130-082-50*

Light protection cover size 3

Protective cover for image plate, format 2.7 x 5.4 cm, 100 pcs

2130-083-00

Light protection cover size 4

Protective cover for image plate, format 5.7 x 7.6 cm, 100 pcs

2130-084-00 VistaScan Combi View/Mini View/Mini Plus

2140-074-00 VistaScan Perio Plus

Bite block protector size 4

Bite block protector for image plates size 4 (occlusal), 100 pcs

2130-074-03

*Also available in white

Simple handling



Image Plate and Film Holder System Set Color

Image plate holders from Dürr Dental offer optimum protection for the image plate with their smoothed and rounded edges. The special stepped design of the right-angle technique Color aiming rings enables simple and safe aiming of the X-ray tube. Holders are available for all applications – their use helps to ensure that accurate X-ray images are produced. The colour coding of the various image plate holders enables simple and error-free handling.

Image Plate and Film Holder System Set Color:

Color indicator arm anterior (1 pce)
Color indicator arm posterior (2 pcs)
Color indicator arm bitewing (1 pce)
Color aiming ring anterior (1 pce)
Color aiming ring posterior (2 pcs)
Color aiming ring bitewing (1 pce)
Color bite blocks anterior (8 pcs)
Color bite blocks posterior (8 pcs)
Color bite wing bite blocks 3 x 4 (2 pcs)
Color bitewing bite block 2.7 x 5.4 (1 pce)
Color bitewing bite block 2 x 3 (1 pce)
2130100015

Image Plate and Film Holder System Set Color Endo

Upgrade Set:

Endo Color indicator arm (1 pce)
Endo Color bite block left anterior (1 pce)
Endo Color bite block left posterior (1 pce)
Endo Color bite block right anterior (1 pce)
Endo Color bite block right posterior (1 pce)
Endo Color aiming ring (1 pce)

2130100014

Disinfection & Cleaning

FD 350, FD 333 disinfection wipes

Particularly suitable for wipe disinfection and cleaning of medical devices. They can be used e.g. for covers for intraoral plates and films, straight and contra-angle hand pieces etc. The disinfection wipes FD 333 offer a fully viricidal solution.

FD 350 Disinfection wipes

12x containers with 110 wipes/carton

CDF35CA0140

FD 333 Disinfection wipes

12x containers with 110 wipes/carton

CDF333A0140




Image plate cleaning wipes

VistaScan image plate cleaning wipes provide intensive cleaning of the image plates and also have a disinfecting effect. Advantage: long-lasting, optimum image quality.

Image plate cleaning wipes

10 wipes

CCB351B1001

