

LumiContrast

Low reflection - High contrast
Ideal for Loupes and Microscopes!







A WORLD FIRST!



SWISS
QUALITY

Order information:

Introkit:	Item No.	Qty
LumiContrast Kit - 120 Matrices + 3 Rings + 6 Delta Silicone Tubes for Rings	REF 6700	

Refill packs:	Item No.	Qty
LumiContrast - Sectional Matrices - Premolar 0.04mm/5mm, 35 pcs 	REF 6705	
LumiContrast - Sectional Matrices - Molar 0.04mm/6.4mm, 35 pcs 	REF 6706	
LumiContrast - Sectional Matrices - Premolar 0.025mm/5mm, 35 pcs 	REF 6735	
LumiContrast - Sectional Matrices - Molar 0.025mm/6.4mm, 35 pcs 	REF 6736	
LumiContrast - Sectional Matrices 0.04mm/6.4mm+Ext., 35 pcs 	REF 6746	
LumiContrast - Sectional Matrices 0.025mm/6.4mm+Ext., 35 pcs 	REF 6756	
LumiContrast - 3 Rings + 6 Delta Silicone Tubes for Rings 	REF 6757	

Address and signature:


My dealer is:

Send your order by fax to:
+41 91 946 32 03

or by e-mail to:
sales@polydentia.com

LumiContrast

A WORLD FIRST!



With usual sectional matrix



With LumiContrast



Tired eyes? LumiContrast is the solution!

Forget your tiring sectional matrices, step into LumiContrast age!

LumiContrast is an **innovative sectional matrix**, result of years of research and development on materials for high standards in dentistry, together with dental specialists from all over the world.

We have considered parameters as **light („Lumi“)** and **contrast („Contrast“)** to develop the best solution for your eyes: an innovative sectional matrix, which **protects your eyes** from unnecessary stress caused by reflection and lack of contrast.

Features and benefits:

- Mat, blue matrix surface for **low reflection** and **no blinding effect**,
- **High contrast** between dark colour/white teeth for precise identification of working area,
- High-Tec matrix colouration method; colour **cannot be scratched** away with handling;
- No colour-coating: **no particles** scattered into the filling material or cavity,
- Premium steel quality, **doesn't collapse into the cavity**, no deformation by wedges.

LumiContrast are a must-have when working with modern technologies with intense light sources, such as loupes and dental microscopes.