

Caps Dispenser

Direct application dispenser for dispensing restorative materials in Caps



REF 9301 Dispenser – Caps

Dispenser type 2

Dispenser to extrude dental materials from cartridges



For 1:1 or 1:2 cartridges (50 ml)

REF 2158 Dispenser – type 2

Dispenser type 3

Dispenser to extrude dental materials from cartridges



For 1:1 or 1:2 cartridges (25 ml)

REF 2259 Dispenser – type 3

QuickMix Dispenser

Dispenser to extrude dental materials from QuickMix syringes



For 1:1 QuickMix syringes

REF 2261 Dispenser – QuickMix

AC Applicator type 1

Applicator for the direct administration of glass ionomer cement restorative and luting materials in application capsules (AC)



REF 2331 Applicator – AC type 1

AC Applicator type 2

Applicator for the direct administration of glass ionomer cement restorative and luting materials in application capsules (AC)



REF 2334 Applicator – AC type 2

AC Activator

Device to activate VOCO application capsules (AC)



REF 9300 Activator – AC

Application brush dispenser

Dispenser for hygienic removal of application brushes – can be repeatedly refilled with two types of brushes



REF 2244 Dispenser – application brushes