



## Lower Molar Clamps



### Rubber Dam Clamp #24

Flat jaws for upper right and lower left molars; paired with RDCM25.

Black Line **RDCM24X**

Satin Steel™ Alloy **RDCM24**



### Rubber Dam Clamp #25

Flat jaws for upper left and lower right molars; paired with RDCM24.

Black Line **RDCM25X**

Satin Steel™ Alloy **RDCM25**



### Rubber Dam Clamp #3

Flat narrow jaws for small lower molars.

Black Line **RDCM3X**

Satin Steel™ Alloy **RDCM3**



### Rubber Dam Clamp #W3

Flat narrow jaws for small lower molars.

Black Line **RDCMW3X**

Satin Steel™ Alloy **RDCMW3**




### Rubber Dam Clamp #7

Flat jaws for lower molars.

Black Line **RDCM7X**

Satin Steel™ Alloy **RDCM7**



### Rubber Dam Clamp #7A

Slightly reversed festooned jaws for large lower molars.

Black Line **RDCM7AX**

Satin Steel™ Alloy **RDCM7A**



### Rubber Dam Clamp #W7

Flat jaws for lower molars.

Black Line **RDCMW7X**

Satin Steel™ Alloy **RDCMW7**



### Rubber Dam Clamp #12A

Serrated offset jaws for lower right molars; paired with RDCM13A.

Black Line **RDCM12AX**

Satin Steel™ Alloy **RDCM12A**