

REBASE II FAST

The classic – tried and tested for years and recognised for hard relinings.



Time and cost saving

The material can also be used chair-side. No waiting time for the patient, no laboratory costs for the dentist.


Pleasant

Low temperature development, no unpleasant taste or smell and MMA-free.

Permanent

A special hardener removes the top oxygen inhibition layer and thus causes complete curing. No breeding ground for bacteria or bad odours – for long-lasting use.

Find out more:



Keyfacts relining

- > Can be used chairside or in the laboratory
- > Hard or soft relining
- > Long durability and reliability