



## Kwik-Bite™

### X-ray holders for horizontal bite wings

The Kwik-Bite™ is a film and phosphor plate holder for horizontal bite wings (standard bite wings). For assessing approximal caries, overhangs and gaps, crown margins and calculus in the interdental spaces, and also the state of the alveolar ridge.

#### Compatibility

The Kwik-Bite is compatible with the most commonly used X-ray film formats (1, 2, 3). The Kwik-Bite with Ring is also suitable for indirect digital systems. It is compatible with all major brands of phosphor plate sensors.

#### Kwik-Bite with Ring


REF/Part No. 1780

5 Kwik-Bite with Ring, 1 Centring Device

#### Kwik-Bite with Index

REF/Part No. 270

15 Kwik-Bite with index



Kwik-Bite with Index



Kwik-Bite with Ring

## Super-Bite™

### X-ray holders for periapical radiographs

The Super-Bite™ is a film and phosphor plate holder indicated for periapical radiographs of anterior and posterior teeth. The system includes two types of holder – anterior (green) and posterior (red) – for use in all four quadrants.

#### Compatibility

The Super-Bite with Ring is also suitable for indirect digital systems. It is compatible with all major brands of phosphor plate sensors.

#### Super-Bite with Ring Assorted Kit

REF/Part No. 1020

2 Super-Bite Anterior, 2 Super-Bite Posterior, 2 Centring Devices

#### Super-Bite with Ring Refill Packs

REF/Part No.

4 Super-Bite Anterior with Ring, 1 Centring Device

1021

4 Super-Bite Posterior with Ring, 1 Centring Device

1022

#### Super-Bite with Index Assorted Kit

REF/Part No. 670

4 Super-Bite Anterior, 4 Super-Bite Posterior, 8 Index

#### Super-Bite with Index Refill Packs

REF/Part No.

8 Super-Bite Anterior, 8 Index

671

8 Super-Bite Posterior, 8 Index

672



Super-Bite Anterior



Super-Bite Posterior