


Cylinder flat end



837




Shank	Ring	Length mm		8,0	8,0	8,0	8,0	8,0
		ISO	Product Nr.	Ø mm				
FG	Black	806 314 111 544 ...	SG837. ...	010	012	014	016	018
	Green	805 314 111 534 ...	G837. ...	010	012	014	016	018
	Blue	805 314 111 524 ...	837. ...	010	012	014	016	018
	Red	805 314 111 514 ...	F837. ...	010	012	014	016	018

Cylinder round edge



837KR



Shank	Ring	Length mm		8,0	8,0	8,0	8,0	8,0
		ISO	Product Nr.	Ø mm				
FG	Green	806 314 158 534 ...	G837KR. ...	010	012	014	016	018
	Blue	806 314 158 524 ...	837KR. ...	010	012	014	016	018
	Red	806 314 158 514 ...	F837KR. ...	010	012	014	016	018

837 Cylinder flat end
Cavity preparation

