

# ESTELITE COLOR

Expressive stains for individual characterisation of the respective restoration.



## Colours

9 different shades make teeth look younger or older. Even fissures or pits can be imitated naturally.

## Opaque

4 opaque shades provide reliable coverage of heavily discoloured areas. Also suitable for covering metal or for repairs.

## Flowability

Due to their consistency, the stains are easy to apply and can be individually shaped with a small probe or the TOKUYAMA brushes.

Find out more:



ESTELITE COLOR is applied evenly over the discoloured areas to raise the colour value.




After light curing, a layer of dentin shade with the desired shade and saturation is applied over the opaquer.



After light curing, the enamel shade is applied.

## Farben - 13:



Customers also bought this additional product:  
Brush No.24

More at >> [www.tokuyama-dental.eu](http://www.tokuyama-dental.eu)