



Gracey Micro Mini Five™ Currettes


Gracey Micro Mini Five™ Currettes feature an elongated shank to access deep, narrow periodontal pockets and a blade that is 50% shorter and 28% thinner than the Standard Gracey Curette.




Gracey 1/2 Micro Mini Five™ Curette	
Harmony	SMS1/2XE2
#9	SMS1/29E2



Gracey 7/8 Micro Mini Five™ Curette	
Harmony	SMS7/8XE2
#9	SMS7/89E2



Gracey 11/12 Micro Mini Five™ Curette	
Harmony	SMS11/12XE2
#9	SMS11/129E2



Gracey 13/14 Micro Mini Five™ Curette	
Harmony	SMS13/14XE2
#9	SMS13/149E2

The Clayton Gracey Story

In the late 1940s, an innovative dentist from the University of Michigan had a vision. Clayton Gracey wanted to “give every dentist the possibility to treat even the deepest and least accessible periodontal pockets simply and without traumatic stretching of the gingiva.”

Dr. Gracey brought his idea to Hugo Friedman, the founder of Hu-Friedy, who supported his concept. Together they developed a series of 14 single-ended area-specific instruments for the successful removal of deposits from root surfaces.

Dr. Gracey’s designs laid the groundwork for the new instruments of today. Currently, Gracey curettes are available as double-ended instruments with multiple handle choices.

The Gracey curettes have come a long way since the original design. Hu-Friedy currently manufactures a variety of styles and patterns, such as After Five® curettes with a longer terminal shank and thinned blade, Mini Five® curettes with After Five® features and a shorter blade, and the 15/16 and 17/18 patterns. An optional Rigid shank is available on all of the Gracey designs.



Hu-Friedy is continually looking for new ways to improve our existing products. Our success is dedicated to the pioneers in dentistry, including Dr. Clayton Gracey, who helped to make a vision become not only a reality, but also a part of every dental office today.