

Hi-Tray Metal

Metal impression trays

Impression systems / Accessories

Zhermack offers a wide range of impression trays for fixed prosthesis, in metal and in plastic, and for removable prosthesis, in rigid plastic. The Zhermack range of impression trays is also suitable for many uses such as models in the practice, bruxism bites, whitening masks, sports masks for bites, impressions for custom impression trays, preliminary impressions, impressions for in-lay and on-lay, etc.

The **Hi-Tray Metal** are BSA stainless steel trays with 16.70 % chrome for standard applications.

Characteristics

- Can be sterilized in a cold solution or in autoclave
- Available in two versions: solid and perforated

Advantages

- High deformation resistance
- Greater precision of the impression
- Rim Lock anchoring system for a perfect adhesion of the material to the impression tray
- Excellent finishing
- Soldering point accuracy
- Durability

Single packaging in the solid impression tray version:

Size	1	2	3	4	5
Upper	D51SS	D52SS	D53SS	D54SS	D55SS
Lower	D51IS	D52IS	D53IS	D54IS	D55IS

Single packaging in the perforated impression tray version:

Size	1	2	3	4	5
Upper	D51SP	D52SP	D53SP	D54SP	D55SP
Lower	D51IP	D52IP	D53IP	D54IP	D55IP

Kit pack: 1 kit contains 5 upper trays and 5 lower trays

D5K10S	solid
D5K10P	perforated

Kit pack: 1 kit contains 4 upper trays and 4 lower trays

D5K08S	solid
D5K08P	perforated



Hi-Tray Metal perforated upper



Hi-Tray Metal solid upper



Hi-Tray Metal perforated lower



Hi-Tray Metal solid lower