



Gracey Curettes


Gracey curettes are area-specific curettes with blades having a single cutting edge, offset at 70° to the shank, to provide a perfect working angulation to the tooth surface.



Learn more about our
[Gracey Curettes](#).

Standard Gracey Curettes

Standard Gracey Curettes allow for deep scaling, root planing and periodontal debridement. Gracey Rigid curettes have thicker, stronger and less flexible shanks than Standard Graceys to allow for removal of moderate to heavy calculus. Gracey Extra Rigid Curettes are ideal for tenacious calculus removal. Standard, Rigid and Extra Rigid Gracey Curettes all have identical blade dimensions.



Gracey 1/2 Curette	Standard	Rigid
Harmony	SG1/2XE2	SG1/2RXE2
#9	SG1/29E2	SG1/2R9E2
#G9, Steel Grip	SG1/2G9E2	SG1/2RG9E2
#C8, Gray	SG1/2C8E2	SG1/2RC8E2
#6	SG1/26	SG1/2R6
#7	SG1/27	SG1/2R7
#8	SG1/28	SG1/2R8
#2	SG1/2	SGR1/2R



Gracey 3/4 Curette	Standard	Rigid
Harmony	SG3/4XE2	SG3/4RXE2
#9	SG3/49E2	SG3/4R9E2
#C8, Gray	SG3/4C8E2	SG3/4RC8E2
#6	SG3/46	SG3/4R6
#7	SG3/47	SG3/4R7
#8, Gray	SG3/48	SG3/4R8
#2	SG3/4	SGR3/4R



Gracey 5/6 Curette	Standard	Rigid	Extra Rigid
Harmony	SG5/6XE2	SG5/6RXE2	N/A
#9	SG5/69E2	SG5/6R9E2	N/A
#G9, Steel Grip	SG5/6G9E2	N/A	N/A
#C8, Yellow	SG5/6C8E2	SG5/6RC8E2	N/A
#6	SG5/66	SG5/6R6	SG5/6ER6
#7	SG5/67	SG5/6R7	SG5/6ER7
#8, Gray	SG5/68	SG5/6R8	SG5/6ER8
#2	SG5/6	SGR5/6R	N/A