

Zimmer® Collagen Plug, Tape and Patch

Protect. Heal. Repair.



ABSORBABLE WOUND DRESSINGS

1 Protects Wound Bed
Adheres and provides coverage to oral wounds and sores

2 Aids Wound Healing
Porous, absorbable matrix supports delicate new tissue

3 Highly Porous,
Highly Absorbent
Zimmer Collagen Plug holds up to 30x own weight in fluid
Resorbs in fewer than 30 days
Greater than 90% open pores

4 Versatile For
Everyday Use
Three convenient shapes for common oral wounds and procedures

ABSORBABLE COLLAGEN MATRICES IN THREE CONVENIENT SHAPES

Zimmer Collagen Plug



10 mm x 20 mm
Cat. No. 0102Z

Zimmer Collagen Tape



2.5 cm x 7.5 cm, 1 mm thick
Cat. No. 0100Z

Zimmer Collagen Patch



2 cm x 4 cm, 3 mm thick
Cat. No. 0101Z

VERSATILE WOUND PROTECTION

- Designed to resorb within a short timeframe
- Indicated for management of oral wounds and sores
 - Denture sores
 - Oral ulcers (non-infected or viral)
 - Periodontal surgical wounds
 - Suture sites
 - Burns
 - Extraction sites
 - Surgical wounds
 - Traumatic wounds



1 Atraumatic tooth extraction.



2 Modified Zimmer Collagen Plug delivered to protect grafted site.



3 Sutured site.



4 Ovate pontic contours soft tissue during healing.

Contact us at 1-800-342-5454 or visit zimmerbiometdental.com

Zimmer Biomet Dental
Global Headquarters
4555 Riverside Drive
Palm Beach Gardens, FL 33410
Phone: +1-561-776-6700
Fax: +1-561-776-1272

Unless otherwise indicated, as referenced herein, all trademarks are the property of Zimmer Biomet; and all products are manufactured by one or more of the dental subsidiaries of Zimmer Biomet Holdings, Inc. and marketed and distributed by Zimmer Biomet Dental and its authorized marketing partners. Zimmer Collagen Plug, Tape and Patch are manufactured by Collagen Matrix, Inc. For additional product information, please refer to the individual product labeling or instructions for use. Product clearance and availability may be limited to certain countries/regions. This material is intended for clinicians only and does not comprise medical advice or recommendations. Distribution to any other recipient is prohibited. This material may not be copied or reprinted without the express written consent of Zimmer Biomet Dental. ZB0030 REV B 03/20 ©2020 Zimmer Biomet. All rights reserved.

