



Chisels

Our bone chisels are meticulously handcrafted from the highest grade surgical stainless steel and engineered to provide optimal performance when reshaping and contouring bone or reflecting secondary palatal flaps.



Learn more about our **Chisels**.

Periodontal Bone Chisels



36/37 Back-Action DE Chisel

#9	C36/379
#6	C36/376

Best in Practice Fact

The Rhodes Back-Action Chisel is designed for use with a pull stroke. It is ideal for removing bone adjacent to the tooth without causing trauma and is especially useful on the distal of last molars.



1/2 Kramer-Nevins Chisel DE

#6	CKN1/26
----	----------------



TG DE Chisel

#6	CTG6
#2	CTG



Modified TG with Ochsenbein Cutouts

#6	TG-O6
#2	TG-O



Urban Chisel

#6	CURBAN1
----	----------------



Urban Bone Chisel

#6	CURBAN2
----	----------------