



STERILIZATION CONTROL: "HELIX TEST" HELIX TEST

- Helix test evaluates the capacity of penetration of the sterilization agent (steam, VH_2O_2 , EO, formaldehyde) inside hollow loads.
- Air and any other non-condensing gases shall be removed from the sterilization chamber before the sterilizing agent comes into contact with the devices to be processed so that the test result is not adversely affected.
- Helix test is suitable for any sterilization process in a loaded autoclave.
- Suitable to be used with all types of autoclaves.
- Full conformity to standards in force.
- N. 1 Helix with 25, 50, 100, 250 strips.

STERILIZATION CONTROL: "EC.STEAM 98.4" STEAM STERILIZATION INTEGRATOR

- Integrators for the control of steam sterilization. Integrators evaluate the effectiveness of the process and exposure to sterilization parameters (time - temperature - saturated steam concentration) to fulfill the steam sterilization cycle at 121°C and 134°C.
- Integrators are completely laminated to prevent ink from migrating to the product and to preserve test result.
- Reliable and easy to use. Non-toxic chemical indicator.
- Full conformity to standards in force.
- Class 4.
- 250 pcs/box.

STERILIZATION CONTROL: "EC.STEAM 98.3-L" STEAM STERILIZATION INTEGRATOR

- Integrators for the control of steam sterilization. Integrators evaluate the effectiveness of the process and exposure to sterilization parameters (time - temperature - saturated steam concentration) to fulfill the steam sterilization cycle at 121°C and 134°C.
- Integrators are completely laminated to prevent ink from migrating to the product and to preserve test result.
- Reliable and easy to use. Non-toxic chemical indicator.
- Full conformity to standards in force.
- Class 5.
- 250 pcs/box.

STERILIZATION CONTROL: "EC.STEAM 98.3L AD" STEAM STERILIZATION ADHESIVE INTEGRATOR

- EC.STEAM 98.3L AD presents the same features as the integrator EC.STEAM 98.3L but, thanks to its adhesive surface, it can be attached after use to the patient's medical record or to sterilization process record sheet.

HELIX TEST - HELIX TEST SYSTEM

REF.	DESCRIPTION	PACK	MIN LOT
HELIX TEST 25 KIT	Helix kit test 25 strip	1 pcs	18 pcs
HELIX TEST 50 KIT	Helix kit test 50 strip	1 pcs	18 pcs
HELIX TEST 100 KIT	Helix kit test 100 strip	1 pcs	18 pcs
HELIX TEST 250 KIT	Helix kit test 250 strip	1 pcs	18 pcs

REF. HELIX TEST SYSTEM



REF. EC.STEAM 98.4



REF. EC.STEAM 98.3-L

