



# Individo® Lux · Profibase

PROVEN TRAY BLANKS AND BASE PLATES FOR  
PARTIAL, FULL AND IMPLANT PROSTHETICS

## Individo® Lux · Profibase

### A SOLID BASE IN THE LABORATORY

#### Individo® Lux

##### Indications

Fabrication of individual trays for use in crown, bridge and total prosthetic techniques

Fabrication of bases for bite patterns

Fabrication of odontograms

##### Advantages

- In blue-transparent or blue-opaque colour
- Allows precise and economic processing
- Stable trays
- Pleasant mint flavour
- Non-sticky consistency
- Light-cured with standard halogen light or UVA lights



##### Presentation

REF 2402	Maxillae 50 pcs. blue opaque
REF 2403	Mandibles 50 pcs. blue opaque
REF 2408	Maxillae 50 pcs. blue transparent
REF 2409	Mandibles 50 pcs. blue transparent

#### Profibase

##### Indications

Fabrication of individual trays for use in crown, bridge and total prosthetic techniques

Fabrication of bases for bite patterns

Base plates for wax contouring in the total prosthetic technique

##### Advantages

- Gingiva shade
- Simple and time-saving shaping by hand
- Permanent stability and compressive strength
- Pleasant mint flavour
- Non-sticky consistency
- Light-cured with standard halogen light or UVA lights



##### Presentation

REF 2111	Maxillae / mandibles 2 × 25 pcs.
REF 2112	Maxillae 50 pcs.
REF 2113	Mandibles 50 pcs.

VOCO GmbH  
Anton-Flettner-Straße 1-3  
27472 Cuxhaven  
Germany  
  
Tel.: +49 (0) 4721-719-0  
Fax: +49 (0) 4721-719-140  
  
info@voco.com  
www.voco.com

Available from: