


**Needle Cap Holders**


Needle cap holds attach to the cassette to help maintain the cap in upright position for one-handed needle recapping. Ask your rep for options to permanently install needle cap holders to your cassettes.

2 Pack **IM1001**

6 Pack **IM1003**

Combination Pack **IM1004**

Includes: 6 snap-in needle cap holders and 6 snap-in hinged instrument clips

**Rubber Dam Clamp Organizing Board**

Rubber Dam Clamp Organizing Board. Holds up to 8 clamps during reprocessing. Fits into cassette accessory of most IMS™ Cassettes.

4.6" X 2.4" X .375" (117 mm x 61 mm x 10 mm) **RDCOB**

Bur Cushions

IMS™ Bur Cushions feature a hollow, rectangular, silicone base with a stainless steel lid designed to hold burs securely in place during cleaning and sterilization. Hollow design allows drainage and proper cleaning. Silicone base allows for different size burs. Please check bur manufacturer instructions for use to ensure ability to be reprocessed.

**IMS™ 1/2 Bur Cushion, Short Lid**

1.75" X 1.3" X .9" (44 mm x 33 mm x 22 mm) **IMS-1372SH**

IMS™ 1/2 Bur Cushion, Tall Lid

1.75" X 2.44" X .9" (44 mm x 62 mm x 22 mm) **IMS-1372TH**

IMS™ Bur Cushion, Short Lid

3.3" X 1.3" X .9" (84 mm x 33 mm x 22 mm) **IMS-1372S**

IMS™ Bur Cushion, Tall Lid

3.3" X 2.44" X .9" (84 mm x 62 mm x 22 mm) **IMS-1372T**